



INTERVENTIONAL  
VASCULAR THERAPY

# Celsite<sup>®</sup> Safety

VENOUS ACCESS PORT SYSTEM

# Celsite® Safety

Access ports for mid to long-term venous applications

Celsite® Safety is intended to be used in any condition that requires mid to long-term intermittent or continuous central venous infusions.

The anatomic design with a low profiled nose simplifies the insertion and allows the creation of a small port pocket to downsize the trauma.


Complete range of lasermarked silicone and PUR catheters with an atraumatic tip and graduated from 5 cm

Solid radiopaque connection ring with anti-kink protection

2 extra large suture holes available with and without silicone plugs to facilitate fixation of the port

Extra large septum diameter for easy puncturing

## PEEK AND TITANIUM COMBINATION AS SPECIAL SAFETY FEATURES



# Celsite® Safety

Long Term Central Venous Access without Compromise on Safety

## LARGE PUNCTURE AREA WITH HIGH DENSITY SILICONE SEPTUM

For simple puncturing and reliable sealing to allow good port life

## EXTRA LARGE SUTURE HOLES

For easy fixation of the access port with sutures

## LASERMARKED PUR AND SILICONE CATHETERS

- Clear readability of the catheter length
- No ink is added

## RADIOPAQUE CONNECTION RING

Anti kink protection and additional fixation of the catheter




## HIGH PRESSURE RESISTANT

- Complete range of Celsite® Safety is resistant to high pressure injections up to 325 psi
- Enables power injections of contrast media
- No need for additional venous access

## RADIOPAQUE CT – MARKING

Clear identification of high pressure resistance under the x-ray



Power  
Injections



325 psi

# Celsite® Safety

Long Term Central Venous Access without Compromise on Safety

## REDUCED TITANIUM CONTENT

To reduce metal induced artifacts in MRI

## PEEK AS HOUSING AND CHAMBER MATERIAL

- Poly Ether Ether Ketone - A biocompatible material with high chemical and pressure resistance and excellent durability characteristics.
- Natural color without any additional substances

## TITANIUM BOTTOM PLATE

High puncture resistance

## TITANIUM EXIT CANNULA

Precision in catheter connection



## Surecan® SAFETY II - PORT NEEDLE

- Intuitive safety mechanism to reduce the risk of needle stick injuries
- High pressure resistance up to 325 psi



## Safecan™ SAFETY - PUNCTURE NEEDLE

- Safety mechanism to reduce the risk of needle stick injuries
- Echogenic puncture needle for needle tip location via ultrasound




# Celsite® Safety

## Accessories

### ACCESSORIES

- A** Tunneling rod for easy catheter tunneling
- B** Peelable introducer sheath with dilatator for easy percutaneous access
- C** Omnifix® Luer Syringe
- D** 1 Surecan® Straight for flushing, aspiration and local anesthesia
- E** Long vein lifter allows easy handling
- F** Separate Rinsing Hub for more flushing flexibility during implantation
- G** Safecan™ Safety - Safety Echogenic Vein Puncture Needle
- H** Surecan® Safety II - Safety Port Needle
- I** Stainless Steel guidewire with flexible J-tip
- J** Ultrasound Cover



	Implantation technique	Seldinger	Surgical cut-down	Seldinger (SNT)
	Kit Designation	1	2	3
A	Tunnelling Rod	✓	-	✓
B	Tear-away Introducer	L 180/140 mm	-	L 180/140 mm
C	Omnifix® Luer Syringe	10 mL	-	10 mL
D	Straight Surecan®	22 G x 30 mm	22G x 30 mm	22 G x 30 mm
E	Vein Lifter	✓	✓	✓
F	Rinsing Hub	✓	✓	✓
G	Safecan™ Safety - Puncture Needle	18 G x 70 mm	-	18 G x 70 mm
H	Surecan® Safety II	20 G x 20 mm	-	-
I	J Guide Wire with Dispenser	0.035" x 50 cm	-	0.035" x 50 cm
J	Ultrasound Cover	-	-	1 US Probe (127x1473,2mm), 1 Gel, 2 Orange Fixation Rings, Pouch

# Celsite® Safety

Celsite® Safety offers a wide range of Silicone and PUR catheters as well as two different port sizes, Standard and Small.

Catheter	OD (F/mm)	ID (mm)	Length (mm)	Flow rate* (ml/min)		Recommended maximum flow rates (mL/s) Contrast media at 37°C (325 psi = 22.4 bar)**			Implantation technique	Type	Reference	Accessories see page 5
						Viscosity up to 11.4 mPa.s (cP)						
						19 G	22 G	19 G				
<b>Standard</b>												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Surgical cut-down	T601F	4437556	2
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Surgical cut-down	T601L	4437573	2
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Surgical cut-down	T601P	4437565	2
PUR	8.5 / 2.8	1.6	500	48	12	2	5	5	Surgical cut-down	T601H	4437581	2
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SST601F	4437603	1
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SST601L	4437612	1
Silicone	10 / 3.2	1.6	500	48	12	2	5	5	Seldinger	SST601G	4437620	1
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Seldinger	SST601P	4437607	1
PUR	8.5 / 2.8	1.6	500	48	12	2	5	5	Seldinger	SST601H	4437617	1
<b>Small</b>												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Surgical cut-down	T605F	4437758	2
Silicone	10 / 3.2	1.6	500	48	12	2	5	5	Surgical cut-down	T605G	4437786	2
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SST605F	4437803	1
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SST605L	4437817	1
Silicone	10 / 3.2	1.6	500	48	12	2	5	5	Seldinger	SST605G	4437822	1
PUR	5 / 1.6	1.1	500	26	10	2	5	5	Seldinger	SST605C	4437800	1
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Seldinger	SST605P	4437809	1
PUR	8.5 / 2.8	1.6	500	48	12	2	5	5	Seldinger	SST605H	4437813	1

\* Gravity infusion of saline (0.9%) through a 22G respectively 19G needle from a height difference of 1 m and a catheter length of 40 cm. According to ISO 10555-1

\*\* Flow rates determined according to ISO 10555-6 with a catheter of 20 cm and Surecan® Safety II 20 mm needle

ALL Celsite® SAFETY PORTS ARE PVC, LATEX AND DEHP FREE

**PVC  
FREE**

**LATEX  
FREE**

**DEHP  
FREE**

Standard	Small
<b>Material:</b> Titanium   PEEK <b>Weight:</b> 8g <b>Internal Volume:</b> 0.5 mL	<b>Material:</b> Titanium   PEEK <b>Weight:</b> 5g <b>Internal Volume:</b> 0.3 mL

# Celsite® Safety

## Additional References with Ultrasound Cover and Silicone Plugs

### Celsite® SAFETY SNT WITH ULTRASOUND COVER\*

Catheter	OD (F/mm)	ID (mm)	Length (mm)	Flow rate** (ml/min)		Recommended maximum flow rates (mL/s) Contrast media at 37°C (325 psi = 22.4 bar)***			Implantation technique	Type	Reference	Accessories see page 5
				19 G	22 G	Viscosity up to 11.4 mPa.s (cP)						
						22 G	20 G	19 G				
<b>Standard</b>												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SNT601F (US Probe)	4437592	3
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SNT601L (US Probe)	4437593	3
<b>Small</b>												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SNT605F (US Probe)	4437594	3
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SNT605L (US Probe)	4437595	3

\* Available in CE marked countries

\*\* Gravity infusion of saline (0.9%) through a 22G respectively 19G needle from a height difference of 1 m and a catheter length of 40 cm. According to ISO 10555-1

\*\*\* Flow rates determined according to ISO 10555-6 with a catheter of 20 cm and Surecan® Safety II 20 mm needle

### Celsite® SAFETY WITH SILICONE PLUGS

Catheter	OD (F/mm)	ID (mm)	Length (mm)	Flow rate* (ml/min)		Recommended maximum flow rates (mL/s) Contrast media at 37°C (325 psi = 22.4 bar)**			Implantation technique	Type	Reference	Accessories see page 5
				19 G	22 G	Viscosity up to 11.4 mPa.s (cP)						
						22 G	20 G	19 G				
<b>Standard</b>												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SST701F	4437605	1
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SST701L	4437614	1
Silicone	10 / 3.2	1.6	500	48	12	2	5	5	Seldinger	SST701G	4437621	1
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Seldinger	SST701P	4437609	1
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Surgical cut-down	T701F	4437560	2
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Surgical cut-down	T701L	4437578	2
<b>Small</b>												
Silicone	6.5 / 2.2	1.1	500	26	10	2	5	5	Seldinger	SST705F	4437805	1
Silicone	8.5 / 2.8	1.2	500	34	11	2	5	5	Seldinger	SST705L	4437818	1
Silicone	10 / 3.2	1.6	500	48	12	2	5	5	Seldinger	SST705G	4437790	1
PUR	5 / 1.6	1.1	500	26	10	2	5	5	Seldinger	SST705C	4437801	1
PUR	6.5 / 2.1	1.4	500	37	12	2	5	5	Seldinger	SST705P	4437807	1
PUR	8.5 / 2.8	1.6	500	48	12	2	5	5	Seldinger	SST705H	4437815	1

\* Gravity infusion of saline (0.9%) through a 22G respectively 19G needle from a height difference of 1 m and a catheter length of 40 cm. According to ISO 10555-1

\*\* Flow rates determined according to ISO 10555-6 with a catheter of 20 cm and Surecan® Safety II 20 mm needle

**Distributor**

**B. Braun Melsungen AG | Vascular Systems | Sieversufer 8 | 12359 Berlin | Germany**  
**Phone +49 30 568207-300 | Fax +49 30 568207-130 | [www.bbraun.com](http://www.bbraun.com)**

Manufacturer acc. to MDD 93/42/EEC: B. Braun Médical, 26 rue Armengaud, 92210 Saint-Cloud, France, [www.bbraun.fr](http://www.bbraun.fr)

The product trademarks „Celsite“, „Omnifix“ and „Surecan“ are registered trademarks of B. Braun Melsungen AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.